
Maggot Debridement Therapy (MDT)

Effective and Accessible Bio Therapeutic Solutions

Company Profile

Cuprina Pte Ltd is a Singapore-based MedTech company. Our first products to market were bio-dressings made of live, sterile maggots.

We launched our MEDIFLY product range in 2020 and they are used across hospitals and veterinary clinics in Singapore to treat chronic wounds.

The MEDIFLY products are specifically used to effectively debride and manage infections associated with chronic wounds.



What is Maggot Debridement Therapy?



Maggot Debridement Therapy (MDT) is the use of sterile, medical-grade maggots (fly larvae) to clean and disinfect chronic wounds. MDT can be used for the removal of necrotic tissue, biofilm, or slough and on bacteria-infested wounds which may be antibiotic-resistant (i.e.: MRSA).

How Maggots Work?

- Maggots are living creatures which require oxygen and liquefied food to survive
- Maggots have no teeth... only “sucking” parts,
- As the maggots crawl on the wound bed, the exo-spines of the maggot work like “mini scalpels” which mimics the action of surgical debridement.
- Maggots are chemical factories - secreting a powerful mixture of proteolytic enzymes which break down dead tissue, liquidizing it as they move on the surface of the wound.
- Maggots then ‘suck’ up this liquidized tissue, ingesting it and digesting it.
- Maggot Enzymes only liquefy devitalized tissue including MRSA and leave healthy tissue undisturbed.
- Maggots also ingest and digest bacteria within the devitalized tissue from the wound, this is done in their gut.
- Secrete chemicals with inherent antimicrobial activity which help combat infections.
- Reduce infections decreases inflammation and promotes wound healing.





Treatment of

ANY CHRONIC ULCER: Pressure Sores, Diabetic Foot Ulcers, Venous, Ischemic, Malignant, Burns....


Methicillin-resistant Staphylococcus aureus (MRSA)

Devitalized tissue – slough, necrosis, gangrene

Non-aggressive & quick debridement

Biofilm formation

Painful adhered Slough



Diabetic Foot Ulcers



Pressure Sores



Slough Wound



Gangrene

Superiority On the 3 Aspects

1. Proven Effectiveness (TIME)

2. Safe & Easy

3. Cost-effective



Wound Bed Preparation (WBP)

The concept of wound bed preparation (WBP) as a clinical approach to the treatment of chronic wounds has been adopted and established as an integral part of the treatment of chronic wounds.

WBP is a holistic concept, which addresses the factors that contribute to the chronicity of wounds, including underlying diseases, and identifies measures to remove barriers to healing.

These measures aim to remove dead tissue and slough, normalise inflammation, re-establish moisture balance and support the movement and migration of cells essential for wound repair.



Effectiveness of MDT (TIME Principles)

T

Tissue: non-viable or deficient

1. Debridement and concomitant loss of bacterial bioburden;
2. Removal of biofilm;
3. Removal of tissue docking sites for bacteria, and effects on bacterial adhesins.

I

Infection or Inflammation

1. Raised wound pH;
2. Antibacterial and antifungal effects;
3. Inhibition of the complement system, and inflammatory cell migration and activation

M

Moisture imbalance

1. larval feeding is facilitated by the liquefaction of tissues resulting from the action of digestive enzymes, which in turn can contribute to the moisture profile of the wound.

E

Edge of wound: non-advancing or undermining

1. Promotion of cell motility and angiogenesis

Tissue - Targeted Debridement

Maggots perform debridement by physically feeding on dead tissue, cellular debris and exudate present in sloughy wounds. Their feeding action physically breaks up necrotic or sloughy tissue, which is then consumed and digested.

This process is mediated by proteolytic enzymatic digestion through a process of extracorporeal digestion. Collagenases, trypsin-like and chymotrypsin-like enzymes are secreted which breakdown devitalised tissue into a semiliquid form which the larvae can ingest.

The larvae of *Lucilia sericata* do not digest living human tissue. This selective process is one of the major advantages of larval debridement therapy as it spares the healthy tissues necessary for healing (Gottrup and Jørgensen, 2011; Fforwm and Meinwe, 2013).




Figure 1. The mean percentage slough pre, during and post MDT (n=11). Error bars indicate standard deviation. There was 45% reduction in slough following the 2nd cycle of MDT compared with baseline ($p < 0.001$).

Infection - Antimicrobial Effect

Additional mechanisms of action include an antibacterial effect within the wound as bacteria contained in liquefied material is ingested and digested, reducing the bioburden within the wound and larval secretions that prevent the formation of, and reduce pre-formed biofilms (Harris et al, 2007; Cazander et al, 2009).

One clinical trial reported that the number of infected wounds decreased between days 1 and 15 with MDT, but not conventional treatment (Mudge et al., 2014). Infections cure faster, and patients remain infection-free for a longer period with MDT (Shi & Shofler, 2014)




Figure 1*. Antibiotic-free days in the maggot debridement therapy and control groups. Horizontal lines represent mean; error bars, SEM; shaded boxes, first to third quartiles. *P = .0001.

Moisture - Rebalance of Humidity

As maggot feeding is facilitated by the liquefaction of tissues resulting from the action of digestive enzymes, this in turn can contribute to the moisture profile of the wound.

At the initiation of MDT, the production of exudate is often enhanced, which is probably a result of the degradation of necrotic and infected material.

As maggots also require a balance of moisture and humidity, they have the potential to act as visual bio-indicators of the wound environment, as they can drown or desiccate in conditions adverse to their survival.



Epithelial Edge - Healing Progression

Two studies [1,2] suggested that MDT hastens the growth of granulation tissue. Sherman reported that ulcers treated with MDT attained at least 50% of granulation tissue within three weeks, whereas Sherman found that the amount of healthy granulation tissue was statistically significant after complete debridement was achieved in four weeks.

Ulcer healing rates are seven times higher with MDT than with conventional therapy (Wilasrusmee et al., 2014)

MDT achieved significantly better granulation rates and shorter healing times than traditional dressings in pressure ulcers and diabetic foot ulcers (Shi & Shofler, 2014)




FIGURE 3. Average percentage of wound base covered by granulation tissue during MDT ($N = 43$) or conventional therapy only ($N = 49$). Error bars indicate standard error; asterisks indicate significant differences in mean percentage of granulation tissue ($p < 0.05$).

Amputation Prevention

Maggot Debridement Therapy is associated with a 3 x lower rate of amputation.[1]

Studies show that when compared to standard therapies. Maggot Debridement Therapy was not only associated with lower amputation rates but also increased healing rates and in some cases patients avoided imminent amputation. [2, 3]

Another case report done in Singapore showed that the limb salvage rate was 90.9% (10/11), while the wound closure rate was 45.5% (5/11) in this challenging patient population. [4]




Figure 1. Comparison of amputation rates [3]

Source: [1] Armstrong et al. 2005. Maggot therapy in "Lower Extremity Hospice" Wound care. Journal of the American Podiatric Medical Association 95 (3), 254 - 257

[2] Gilead L, Mumcuoglu KY, Ingber A. 2012. The use of maggot debridement therapy in the treatment of chronic wounds in hospitalized and ambulatory patients. J Wound Care; 21(2):78-85.

[3] Tian X et al. 2013. Maggot debridement therapy for the treatment of diabetic foot ulcers: a metaanalysis. Journal of Wound care; 22(9): 462-9.

[4] Reintroduction of maggot debridement therapy in the treatment of diabetic foot ulcers in Singapore: a single institution's initial experience

Safe & Easy

MDT is relatively safer and has lower side effects. The results [1] were promising with minimal side effects besides discomfort and mean visual analogue pain score was 3.3.

Food & Drug Administration (FDA) of the United States has confirmed the prescription and use of Maggot therapy and enacted regulations for it. Also, in other countries, MDT has been used as a kind of drug and at least 24 laboratories in more than 30 countries in the world in 2009 have worked to prepare therapeutic maggots [2].

The MDT application of the takes the same length of time as a standard wound care dressing and can be removed from the wound quickly and easily [3].



Mr William Teo (left) underwent five cycles of maggot debridement therapy at the National University Hospital. Photo: Ili Nadhirah Mansor/TODAY

Source: [1] Wounds Asia 2022 | Vol 5 Issue 2 | ©Wounds Asia 2022 | www.woundsasia.com

[2] Sherman RA. Maggot therapy takes us back to the future of wound care: new and improved maggot therapy for the 21st century. *Journal of diabetes science and technology* 2009;3(2):336-44.

[3] Bültemann. 2015. Anleitung zur Applikation medizinischer Fliegenlarven (*Lucilia sericata*) bei der Behandlung chronischer Wunden. *Wund Management* 2015; 9 (3), S. 100-106



Cost-effective

The median cost of Maggot Debridement Therapy (6,700 Baht) was even less than half of the expenditure for conventional therapy (16,800 Baht). [1]

The cost saving by the use of MDT was due to its shorter curing duration, about half of the curing duration for conventional therapy. [1]


The study has demonstrated that based on the available evidence MDT is shown to be less costly and more effective than the other debridement methods tested. [2]

Costs	Maggot therapy cost (Thai Baht), (n = 80)	Conventional therapy cost (Thai Baht), (n = 70)
Main therapeutic cost		
- Median cost	3,000	8,400
- Range	6,750	14,400
(Min. – Max.)	(2,250 – 9,000)	(3,000 – 17,400)
- Interquartile range	3,000 – 6,000	4,200 – 12,000
Related cost		
- Median cost	2,960	8,400
- Range	15,330	17,400
(Min. – Max.)	(210 – 15,540)	(0 – 17,400)
- Interquartile range	210 – 8,140	4,100 – 11,100
Total cost		
- Median cost	6,700	16,800
- Range	16,080	30,300
(Min. – Max.)	(2,460 – 18,540)	(4,500 – 34,800)
- Interquartile range	3,210 – 14,140	8,400 – 22,200

Table 1. Cost comparison between maggot therapy and conventional therapy

Curing duration	Maggot therapy (Days), (n = 80)	Conventional therapy (Days), (n = 70)
- Median curing duration	6.5	14
- Range	18(3-21)	40(5-45)
(Min. – Max.)		
- Interquartile range	3-11	7-20

Table 2: Comparison of curing duration in maggot therapy and conventional therapy



Source: [1] Eamkong, Suwannee & Pongpanich, Sathirakorn & Rojanaworarit, Chanapong. (2010). COMPARISON OF CURING COSTS BETWEEN MAGGOT AND CONVENTIONAL THERAPIES FOR CHRONIC WOUND CARE. *Journal of Health Research*.

[2] Bennett H. et al 2013. Cost-Effectiveness of Interventions for Chronic Wound Debridement; an Evaluation in Search of Data. *London: Wounds UK 2013; 9(4) Suppl.*



Efficacy

- Targeted Debridement
- Promote Growth of Granulation Tissue
- Antimicrobial Effect
- Reduce Amputation Rate
- Effective for Complex and Hard-to-reach Wounds

Safety

- Minimal Side Effects
- Low Pain Levels
- Widely Used in the World
- Approved by Health Sciences Authority (HSA) Singapore

Price

- Cost-effective
- Shorten Length of Hospital Stay

Our MDT Products



Free Range

- Applied for maximum of 72-hrs
- Allows natural mechanical benefit of maggots' movement
- Extensive coverage



Baggots

- Applied for maximum of 96-hrs
- Reduces pain
- Removes 'Yuck' factor
- Time management – quick application and removal



MDT Application



Step 1

Primary Dressing

Frame the wound with Hydrocolloid dressing



Step 2

Primary Dressing

Place live maggots onto gauze and invert onto the wound



Step 3

Primary Dressing

Encage using gas/air permeable tape e.g. Tegaderm



Step 4

Secondary Dressing

Place gauze slightly above the Bio Dressing

Case Study (1/3)

- 68 years old female patient.
- Post surgical infection.
- Methicillin-resistant Staphylococcus aureus, MRSA.
- Recommended forefoot amputation.



Day 1: Initial presentation chronic MRSA wound



Day 3: Post 1st application of MDT. Necrotic edges are reduced and evidence of granulation present



Day 5: Post 2nd application of MDT. Granulation tissue clearly evident



Day 7: Post 3rd application of MDT. No further vials needed due to nicely granulating base



Day 14: Wound after MDT before discharging

Case Study (2/3)

- 80-year-old diabetic male complicated with nephropathy and neuropathy.
- Clinical history of peripheral arterial disease with mild chronic venous insufficiency.
- Presented with a chronic non-healing ulcer on the left malleolus.
- Sloughy, erythematous, and tender wound.
- Treated with four cycles of Baggot MDT with no adverse effects.



Pre-MDT:
NVT: 100%

1st Removal Change:
NVT: 70%

2nd Removal Change:
NVT: 55%

3rd Removal Change:
NVT: 10%

Post-MDT (Completion)

NVT: Non-Viable Tissue



Case Study (3/3)

- A 59-year-old male presented with worsening gangrene of left posterior heel despite previous revascularization and DAPT (Aspirin + Clopidogrel).
- Methicillin-resistant Staphylococcus aureus, MRSA.
- Wound care management: Alternating treatment involving 1-week MDT, 1- week NPWT and back to 1-week of MDT.



Pre-MDT



3rd Cycle of MDT
&
Start Negative
Pressure Wound
Therapy (NPWT)



5 Days of NPWT
&
Restart MDT



6th Cycle of MDT
&
Restart NPWT



3-Weeks of NPWT



Completion

Frequently Asked Questions (FAQs)

- **Can the maggots turn into flies in my wound?**

No. The next stage of the natural life cycle of a maggot is to form into a pupa before becoming a fly. The maggots can only pupate in dry conditions- not a moist wound.

- **Is it painful?**

Most of the patients on MDT mentioned that they feel no pain or just a very slight pain from the treatment.

- **Are there any side effects?**

No. There are no side effects of MDT.

- **What do care providers for the patient need to do?**

Leave the primary dressing alone until the removal of the maggots. Only change the secondary dressings once every 3-4 hours or whenever the secondary dressing is soaked by the exudate. This prevents the maggots from drowning.

- **How long does the maggot debridement therapy take and how long to leave it for?**

The application of MDT takes 20-25 mins by the nurses or podiatrist. The maggots are left in the wound for 48 hours for the free range and 72 hours for the Baggots(Maggots in a bag).



Attained Regulatory Standards

National Environment Agency
(NEA)



Attained Singapore's NEA License to collect, keep, breed and sell the *Lucilia Cuprina* Blowfly

International Organization for Standardization (ISO 13485)



Quality Management Systems and Manufacturing Processes and Controls have attained ISO 13485 standards

Health Science Authority
(HSA)



Attained Class C Medical Device Classification under the Health Sciences Authority Singapore

Our Customers



Singapore



Malaysia



Hong Kong



瑪麗醫院
Queen Mary Hospital



LKS Faculty of Medicine
The University of Hong Kong
香港大學李嘉誠醫學院



威爾斯親王醫院
PRINCE OF WALES HOSPITAL



東區尤德夫人那打素醫院
PAMELA YUODE NETHERSOLE EASTERN HOSPITAL

Clinical Papers & Media

Case reports

Reintroduction of maggot debridement therapy in the treatment of diabetic foot ulcers in Singapore: a single institution's initial experience

Key words:

- Diabetic foot
- Maggot debridement therapy
- Reintroduction of debridement
- Wound healing

Background: Diabetes mellitus is often associated with peripheral vascular disease and limb ulcers. Among diabetics, the lifetime risk of developing a foot ulcer is approximately 15%. Maggot debridement therapy (MDT) uses sterile, medical-grade maggots for non-surgical treatment of necrotic or sloughy wounds. **Aim:** This study aims to ascertain the safety and efficacy of MDT using locally produced, sterile maggots of *Lucilia sericata* as a means of biological debridement in diabetic foot ulcers (DFUs). **Methods:** We prospectively recruited patients with DFUs to undergo MDT. Baseline and interim wound characteristics, change in slough or necrotic tissue were studied. **Results:** Our results were promising with minimal side effects besides discomfort and mean visual analogue pain score was 3.3. All wounds showed reduction in slough, with an average reduction of 53% after one application and 65% after two applications (p<0.001). **Conclusion:** MDT has proven to be a safe and effective method of debridement, with the pain resulting from MDT easy to manage.

Diabetes mellitus is a condition often associated with peripheral vascular disease and limb ulcers. Diabetic foot ulcers (DFUs) are one of the most prevalent complications in patients who have suboptimal control of their diabetes. There are 1.3 million people with diabetes in Singapore, with 1.1 million being type 2 diabetes. The Singapore Diabetes Prevalence Study (SDPS) reported that 1.1 million Singaporeans are affected by diabetes. In 2010, the Singapore Diabetes Prevalence Study (SDPS) reported that 1.1 million Singaporeans are affected by diabetes. In 2010, the Singapore Diabetes Prevalence Study (SDPS) reported that 1.1 million Singaporeans are affected by diabetes. In 2010, the Singapore Diabetes Prevalence Study (SDPS) reported that 1.1 million Singaporeans are affected by diabetes.

Diabetes mellitus is a condition often associated with peripheral vascular disease and limb ulcers. Diabetic foot ulcers (DFUs) are one of the most prevalent complications in patients who have suboptimal control of their diabetes. There are 1.3 million people with diabetes in Singapore, with 1.1 million being type 2 diabetes. The Singapore Diabetes Prevalence Study (SDPS) reported that 1.1 million Singaporeans are affected by diabetes. In 2010, the Singapore Diabetes Prevalence Study (SDPS) reported that 1.1 million Singaporeans are affected by diabetes. In 2010, the Singapore Diabetes Prevalence Study (SDPS) reported that 1.1 million Singaporeans are affected by diabetes.



Maggot debridement therapy in the tropics – Preliminary outcomes from a tertiary hospital

Nanthakumari Gunasegaran¹, Vivian Qian Hai Seah², Shin Yuh Ang³, Parula Alaworn¹, Wee Ting Gob⁴, Angela Yi Jia Liew⁵, Wei Xian Tan⁶, Helen Yee Ting Tay⁷, The Tee Chong⁸

¹Nanyang Hospital, Singapore General Hospital, Singapore
²Department of Plastic Surgery, Singapore General Hospital, Singapore

ARTICLE INFO

Keywords: Debridement, Maggot therapy, Diabetic foot ulcers, Wound healing, Medical-grade maggot, Pain management, Visual analogue

ABSTRACT

This paper aimed to describe the clinical outcomes and patient acceptance of Maggot debridement therapy (MDT) in a tertiary hospital in Singapore. Patients with non-healing DFUs covering at least 20% of wound bed on three trials with or without debridement or with debridement on the bedside were recruited between January and August 2012. Sterile medical-grade maggots of *Lucilia sericata* were used. Wound episodes were assessed before and after debridement to the type of MDT, either Maggot or Free range MDT, and the number of maggots required to debride the wound. Wounds that were assessed at a tertiary hospital or multiple day patient before the start of MDT therapy, during the wound review or 48 or 72 hours after each cycle of MDT and completion of therapy. Three patients received Maggot therapy, while the remaining 12 received Free therapy. The mean age for patients receiving Maggot MDT was 70.3 (SD = 16) and 63.5 (11.4), respectively. Each patient received three cycles of MDT therapy, on average. The mean common type of wound was open granulating toe wound (4 = 33%), while the least common wound setting was arterial ulceration (2 = 15%). A reduction of 50% was observed in 11 out of 15 patients, and in another 11 patients achieved complete debridement of their DFUs within 10 weeks. Five out of 14 patients had to undergo reoperation within the same admission due to poor wound healing, and three of them were referred back to achieve successful debridement. MDT was well accepted by the patients, and they did not require analgesia. MDT can facilitate wound healing through atraumatic debridement and potentially reduce the need for operative debridement which is deep reprofiling the type of MDT that is optimal to use in tropical countries with high humidity.

1. Introduction

Wound healing is an intricate and regulated process that is crucial in maintaining the protective function of the skin [1]. The restoration of this protective barrier is sometimes required by consolidation [2]. Chronic wounds are the lower limb that cannot be managed by using standard medical grade maggot, may result in the formation of a chronic wound. Features of chronic wounds include persistent inflammation, prolonged inflammation, formation of new blood vessels, and lack of response to repetitive medical [3]. Chronic wounds start as significant barrier on health care, as mentioned and the patient. The Singapore General Hospital (SGH), the main teaching hospital in Singapore, has been using MDT since 2006. This study aims to describe the clinical outcomes and patient acceptance of MDT in a tertiary hospital in Singapore.

2. Methods

This study was conducted in a tertiary hospital in Singapore. Patients with non-healing DFUs covering at least 20% of wound bed on three trials with or without debridement or with debridement on the bedside were recruited between January and August 2012. Sterile medical-grade maggots of *Lucilia sericata* were used. Wound episodes were assessed before and after debridement to the type of MDT, either Maggot or Free range MDT, and the number of maggots required to debride the wound. Wounds that were assessed at a tertiary hospital or multiple day patient before the start of MDT therapy, during the wound review or 48 or 72 hours after each cycle of MDT and completion of therapy. Three patients received Maggot therapy, while the remaining 12 received Free therapy. The mean age for patients receiving Maggot MDT was 70.3 (SD = 16) and 63.5 (11.4), respectively. Each patient received three cycles of MDT therapy, on average. The mean common type of wound was open granulating toe wound (4 = 33%), while the least common wound setting was arterial ulceration (2 = 15%). A reduction of 50% was observed in 11 out of 15 patients, and in another 11 patients achieved complete debridement of their DFUs within 10 weeks. Five out of 14 patients had to undergo reoperation within the same admission due to poor wound healing, and three of them were referred back to achieve successful debridement. MDT was well accepted by the patients, and they did not require analgesia. MDT can facilitate wound healing through atraumatic debridement and potentially reduce the need for operative debridement which is deep reprofiling the type of MDT that is optimal to use in tropical countries with high humidity.

* Corresponding author: Dr. Nanthakumari Gunasegaran, Nanyang Hospital, Nanyang Division, Singapore. Tel: 65 6375 2222. Email address: nanthakumari.gunasegaran@sgg.com.sg (N.G.).
 E-mail address: vivianqianhai.seah@sgg.com.sg (V.Q.H.S.); shinyuh.ang@sgg.com.sg (S.Y.A.); parula.alaworn@sgg.com.sg (P.A.); weeting.gob@sgg.com.sg (W.T.G.); angela.yijia.liew@sgg.com.sg (A.Y.J.L.); weixian.tan@sgg.com.sg (W.X.T.); helenyee.ting.tay@sgg.com.sg (H.Y.T.T.); thetee.chong@sgg.com.sg (T.T.C.).

National University Hospital

The combined use of maggot debridement therapy (MDT) and negative pressure wound therapy for the non-surgical management of diabetes foot osteomyelitis (DFO)

Booy J, Lai A*, Chee YHF, Deng WQ¹
¹Podiatry, Department of Rehabilitation, National University Hospital Singapore
²Department of Orthopaedic, National University Hospital Singapore
³Department of Endocrinology and Nephrology, ChangGung Emergency Medical Centre

CUIPRINA

Introduction

Effective debridement is one of the key principles of wound bed preparation. The presence of bacteria in chronic wounds compromise the effectiveness of regular debridement. Maggot debridement therapy is a form of bio-debridement of non-healing tissue in limbs and acute fire trauma. This research provides evidence suggesting highly successful tissue without damaging healthy tissue.

The "DIG" model of mechanism

Illustrative case

Case 1

A 65-year-old Chinese male presented with worsening gangrene of his plantar foot (Fig. 1A). After debridement, the patient was treated with MDT and NPWT. The patient was discharged after 10 weeks of treatment (Fig. 1B, C). The patient was discharged after 10 weeks of treatment (Fig. 1B, C). The patient was discharged after 10 weeks of treatment (Fig. 1B, C).

Conclusion

Maggot debridement therapy promotes successful healing by removal of necrotic and other cellular debris in patients with chronic wounds. The combined use of MDT and NPWT may be a promising approach for the non-surgical management of diabetes foot osteomyelitis.

More Papers to Read:

- Sherman, R. "Indications, Contraindications, Interactions, and Side-effects of Maggot Therapy." *A Complete Guide to Maggot Therapy: Clinical Practice, Distribution, and Ethics* (2022): 63-78.
- Stadler, F. "The maggot therapy supply chain: a review of the literature and practice." *Medical and veterinary entomology* 34.1 (2020): 1-9.
- Sherman, Ronald A. "Maggot therapy takes us back to the future of wound care: new and improved maggot therapy for the 21st century." *Journal of diabetes science and technology* 3.2 (2009): 336-344.



MEDIA

- CNA - Maggots eating rotting flesh: The unusual hospital treatment that's helping save limbs
- TODAY - Unlikely helpers: 'Fresh' maggots that are hungry to feed on patients' wounds as healing therapy
- 8 Worlds (8视界新闻网) - 晨光善方保健：慢性伤口护理疗法 蛆虫治疗免除截肢

Contact Us



DAVID QUEK
+65 9829 8287
david.quek@cuprina.com.sg

www.cuprina.com.sg

CARL BAPTISTA
+65 8484 4241
carl@cuprina.com.sg